
Thermo Scientific Nicolet iS5N FT-NIR Spectrometer

Near-infrared for today's laboratory

The Thermo Scientific™ Nicolet™ iS™5N FT-NIR spectrometer provides superior performance in a compact size at an affordable price. Combining flexible sample handling and leading Thermo™ Scientific™ OMNIC™ software, the Nicolet iS5N spectrometer sets the benchmark in its class of FT-NIR instrumentation.





Based on the popular Thermo Scientific Nicolet iS5 mid-infrared spectrometer, the Nicolet iS5N FT-NIR spectrometer integrates high-performance optics into a small, rugged package. The spectrometer is perfect for use in a wide range of settings, from small industrial facilities to global manufacturing sites.

The Nicolet iS5N spectrometer withstands hot, humid environments while its user-serviceable components greatly reduce service visits and lower maintenance costs.

The Nicolet iS5N FT-NIR Spectrometer

- Affordable price
- Superior performance
- Small footprint
- Flexible sampling

Unique Optical System

The sealed and desiccated optical unit protects the instrument from humidity and solvent vapors. A self-compensating, dynamically aligned interferometer removes any tilt and shear scanning error, automatically tunes the instrument for best throughput and provides analysis speed for real time survey or screening. Pinned-in-place adjustment-free optics guarantee consistent system performance with minimum maintenance.





Performance

- Spectral range: 11,000–3800 cm⁻¹ optimized, CaF₂ beamsplitter
- **Signal-to-noise:** 10 uABS RMS in less than a minute
- Spectral resolution: Better than 4 cm-1
- Wavelength precision: 0.01 cm⁻¹ at 5300 cm⁻¹
- Linearity: Slope 1 ± 0.05



Easy, user-accessible source replacement

Optics

- **Detector:** Extended range indium gallium arsenide (InGaAs)
- **Source:** High intensity halogen NIR; user replaceable from bottom plate
- Beamsplitter: CaF2 near-infrared
- Laser: Temperature stabilized solid-state Near-IR diode laser

Electronics

 USB 2.0 high-speed bidirectional communication to PC

Dimensions

- **Size:** 35 cm W × 28 cm D × 26 cm H (13.5" W × 10.9" D × 10.2" H)
- Weight: 10 kg (22 lbs)

Regulatory

• Regulatory approvals: CE, ETL

Performance Verification

- Automated performance verification to meet customer ISO/GLP requirements
- External 0.8 mm polystyrene film
- Automatic accessory detection and performance verification
- System Suitability tests for complete system performance assurance

Humidity and Vapor Protection

- Tightly sealed and desiccated optical bench with protective Quartz windows
- Rechargeable desiccant cartridges with humidity indicator
- Internal diagnostics with electronic humidity and heat sensors



High performance optics in a rugged, lightweight package

Diagnostics

- Automatic accessory recognition
- Temperature
- Humidity

Serviceability

- User-replaceable components (without opening cover):
- Source
- Desiccant
- Power supply
- Sample compartment windows

OMNIC Software

OMNIC software for the Nicolet iS5N spectrometer includes these standard features:

- Complete collection of spectral processing and analysis tools
- QCheck correlation for material verification
- Live display of data collection and spectral data preview
- Customizable toolbar and menu options
- User logins and password protection
- Peak analysis tools: peak area, peak height
- Data processing, conversions and corrections
- Spectral Search and Library Manager
- Infrared spectral interpretation with online interpretation guide
- Automatic atmospheric suppression
- · Full-featured report generator
- Quant prediction for PLS, PCR, Beer's Law, CLS, peak ratio, and Discriminant Analysis
- Windows® 7 (32 or 64 bit) and Windows 8 (64 bit) compatible

Minimum PC Requirements

Microsoft® Windows 7 or 8 OS, 2 GB RAM, 3 GHz processor, 160 GB HDD, USB 2.0 or $3.0, 800 \times 600$ CRT/LCD display, DVD drive

www.thermofisher.com/iS5N

©2016 Thermo Fisher Scientific Inc. All rights reserved. Microsoft and Windows are registered trademarks of Microsoft Corporation. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0 Australia +61 3 9757 4300 Austria +43 810 282 206 Belgium +32 53 73 42 41 Canada +1 800 530 8447 China +86 21 6865 4588

Denmark +45 70 23 62 60 Europe-Other +43 1 333 50 34 0 Finland/Norway/Sweden +46 8 556 468 00 France +33 1 60 92 48 00 Germany +49 6103 408 1014 India +91 22 6742 9494 Italy +39 02 950 591 Japan +81 45 453 9100 Latin America +1 561 688 8700 Middle East +43 1 333 50 34 0 Netherlands +31 76 579 55 55



rmo Electron Scientific Instruments LLC,

Russia/CIS +43 1 333 50 34 0 Spain +34 914 845 965 Switzerland +41 61 716 77 00 UK +44 1442 233555 USA +1 800 532 4752



A Thermo Fisher Scientific Brand